



VARIAN

Inspiring Excellence

Background Correction



VARIAN

What is Background Absorption


Non specific absorption caused by:

- Molecular absorption in the Gas Phase
- Light scattering by particles in the Light Path



VARIAN

Molecular Absorption in the UV Region



From data supplied by Dr. John Willis



VARIAN

Background

Radiation from hollow cathode lamp is attenuated by **NON-ATOMIC** source

- Molecular species
- Solid particles
- Absorption
- Scatter

Signal is added to atomic signal

- Results in **FALSELY HIGH SIGNAL**

Most severe in graphite furnace

- Can exceed 2.0 abs



Background / non-specific absorption is normally very small in flame AA

- Chemically rich flame environment dissociates salt particles and molecules very efficiently

What conditions can lead to background absorption in flame AA?

- Analytical wavelengths less than 250 nm
- Low analyte concentrations
- High dissolved solids (salt) solutions
- Fuel rich (cooler) Air/Acetylene flames



VARIAN

General Background Correction

Total absorbance measured

- Atomic + non-specific

Background measured

- Non-specific only

Measurements are time separated

- A few milliseconds

Atomic absorption calculated

- Total absorbance - background absorbance = atomic absorbance

Graphite furnace signals have rapid rise and rapid decay times

- Up to 20 absorbance units/second

Time separated total and background absorbances need to be rapidly made

- Ideally simultaneously
- 2 - 10 ms intervals in commercial instruments
- Larger time difference - greater error
 - Exception is flame AA where we work with a steady state absorption signal

Deuterium
Zeeman
Smith-Hieftje



VARIAN

Deuterium Technique

Most common

Continuum source to measure background

- Deuterium lamp

Operating range from 180 to ~ 425 nm


Background is most significant at shorter wavelength

- Deuterium works well **MOST** of the time



VARIAN


Deuterium Background Correction





VARIAN

Deuterium Background Correction





Radiation from **BOTH** hollow cathode lamp and deuterium lamp are coincident

- If **NOT** measurements made on different atom population
 - Significant error

Hollow cathode lamp energy attenuated by **BOTH** atomic and background species

- Total absorption

Deuterium energy attenuated by background species

- Background only
- Atomic component too small to detect



VARIAN


Calculation of Background

Hollow cathode lamp signal = AA + BGD

Deuterium lamp signal = BGD only

Electronically processed signal = AA only


Single Beam Schematic – D₂





VARIAN

Double Beam Schematic – D₂





Limitations of deuterium background correction

- Intensity of continuum inadequate at high wavelength
- Cannot accurately correct for structured background
- Spectral interferences can occur
 - Rare

Zeeman background correction overcomes these limitations



Atomic spectral lines are split in the presence of a magnetic field


In “simple” or “normal” zeeman effect


- One pi (π) component
 - At original wavelength
 - Polarized parallel to magnetic field
- Two sigma (σ) components
 - Symmetrically displaced around original wavelength
 - Polarized perpendicular to magnetic field



VARIAN

Zeeman Background Correction







VARIAN

Zeeman Background Correction


Total absorbance measured with **magnet OFF**

- Same measurement made by Deuterium or Smith-Hieftje systems

Background absorbance measured with **magnet ON**

- Polarizer excludes pi (π) component
- Measurement made exactly at the analyte wavelength

Molecular species are unaffected by magnetic field



Correction takes place at the **EXACT** analyte wavelength

Correction over the complete wavelength range

Correction for structure background

Correction for some spectral interferences

Only one light source is required

True double beam performance

- Magnet On - magnet Off
- Automatic compensation for lamp drift

Calibration roll-over

- Degree of roll-over is element dependent


Sensitivity loss for some elements

- Degree of sensitivity loss is element dependent
- Expressed as magnetic sensitivity ratio (MSR)
- Majority of elements MSR loss $\leq 10\%$



VARIAN

Zeeman Calibration Roll-Over



Good for difficult samples

- High background
- Unknown interferences

Good when spectral interferences occur

- Se in the presence of high Fe or phosphate
- As in the presence of high Al or phosphate
- Pb in the presence of high phosphate
 - 217.0 nm only