



**VARIAN**


*Inspiring Excellence*

# **AA Light Sources**



VARIAN


# Hollow Cathode Lamp Design





VARIAN


# HCL Operation





VARIAN

# Self-Reversal





VARIAN

# Choosing an Analytical Wavelength

Sample Concentration

Precision

Dilution Factors

Wavelength

- Calibration Curvature

Relative Sensitivity

- Concentration & Dilution

Relative Intensity

- Signal to noise



VARIAN

# Lamp Data

Element	Maximum Current mA	Recommended Current mA	Wavelength	Slit nm	Relative Sensitivity	Relative Intensity
<b>Ni</b>	<b>10</b>	<b>4</b>	<b>232.0</b>	<b>0.2</b>	<b>1</b>	<b>5</b>
			<b>352.5</b>	<b>0.2</b>	<b>5</b>	<b>100</b>
			<b>351.5</b>	<b>0.2</b>	<b>10</b>	<b>30</b>
			<b>362.5</b>	<b>0.2</b>	<b>500</b>	<b>10</b>
<b>Cu</b>	<b>10</b>	<b>4</b>	<b>324.8</b>	<b>0.5</b>	<b>1</b>	<b>100</b>
			<b>327.4</b>	<b>0.5</b>	<b>2</b>	<b>50</b>
			<b>217.9</b>	<b>0.2</b>	<b>8</b>	<b>3</b>
			<b>218.2</b>	<b>0.2</b>	<b>10</b>	<b>2</b>
			<b>222.6</b>	<b>0.2</b>	<b>40</b>	<b>5</b>
			<b>244.2</b>	<b>1</b>	<b>270</b>	<b>15</b>



VARIAN

# Multi Element Lamps

Available for:

Ag-Cd-Pb-Zn

Ag-Cr-Cu-Fe-Ni

Ca-Mg

Ca-Mg-Al

Co-Mo-Pb-Zn

Cu-Fe-Mn-Zn

Cu-Fe-Si-Zn

Cu-Zn

Fe-Co-Ni-Mn-Cu-Cr

Na,K



VARIAN

## Using Multi Element Lamps

Use all normal recommended wavelengths

Refer to Data Sheet for recommended slit width

Refer to Data Sheet for recommended lamp current

The SAME current is used for ALL elements in the lamp – may need to change in the method





VARIAN

# Multi Element Lamps

## Advantages

- More elements can be run without changing lamps
- Lower initial cost

## Disadvantages

- More complex emission profile
- Some secondary lines cannot be used
- Some changes in sensitivity observed
- Operating current may be higher
- Line intensities may be lower

We could increase the lamp current - but this results in the following undesirable changes

- The atom cloud extends beyond excitation region and cools down
- Atoms in the atom cloud may absorb light from the lamp (self-reversal)
- This causes broadening of the emission line and increased calibration curvature


# The High Intensity UltrAA Lamp





VARIAN

# The Varian UltrAA Lamp





2 anodes are fitted inside the UltrAA lamp:

- one for the normal HCL discharge
- the second for a separate boost discharge which is directed across the cathode

Electrons emitted by a heated filament are passed through the atom plume to re-excite the atoms by collision

This increases the emission intensity of the lamp



VARIAN

# UltrAA Lamps - Elements Available

## Single Element UltrAA Lamps

Sb	Cu	Mn	Si
As	Ge	Ni	Te
Bi	Au	Pd	Tl
B	Fe	Pt	Sn
Co	Pb	Se	

## Multi Element UltrAA Lamps

Ag-Cr-Cu-Fe-Ni	Cu-Fe-Mn-Zn
Ag-Cd-Pb-Zn	Cu-Fe-Si-Zn
As-Cu-Fe	Cu-Zn
Ca-Mg-Al	Fe-Co-Ni-Mn-Cu-Cr
Co-Mo-Pb-Zn	Ni-Se


Varian's boosted discharge lamps provide:

- Improvements in intensity
- Extended lifetime
- Improved sensitivity
- Improved calibration linearity



VARIAN

# Lead Lamp Comparison








VARIAN

# Advantages of UltrAA Lamps

Reduced baseline noise

Peak to peak baseline noise for Se is 0.013 Abs. for the conventional HCL compared with 0.008 Abs. for the UltrAA Lamp






VARIAN

# UltrAA Lamp Increased Sensitivity

## Increased sensitivity

- Typical signals for 75  $\mu\text{g/L}$  Se are illustrated






VARIAN

# Advantages of UltrAA Lamps

## Enhanced calibration linearity

- Typical calibrations for Se using standards of 25, 50 and 75  $\mu\text{g/L}$  Se are illustrated



## Improved Characteristic Concentrations and Detection Limits

Element	Wavelength nm	UltrAA Lamp	Conventional HCL
---------	------------------	----------------	---------------------

Characteristic Concentration (in  $\mu\text{g/L}$  with  $20\mu\text{L}$  sample)

As	193.7	0.25 - 0.3	0.31
Pb	283.3	0.15 - 0.26	0.26
Se	196.0	0.49 - 0.5	0.85 - 2.2


Detection Limits (3 sigma in  $\mu\text{g/L}$  with  $20\mu\text{L}$  sample)

As	193.7	0.35 - 0.7	0.3 - 1.4
Pb	283.3	0.18 - 0.23	0.8
Se	196.0	0.28 - 0.57	1.1 - 3.1



VARIAN


# Gold Lamp Comparison





VARIAN


# Comparison of Gold UltrAA and Normal Lamps at Detection Limit Levels





VARIAN

# Lifetime of UltrAA Lamps



## Improved analytical performance

- Lower baseline noise levels
- Increased sensitivity
- Lower detection limits
- Enhanced calibration linearity

Greater reliability

Longer lamp lifetimes